

Carestream

DRX Plus 2530C Detector



# Small. Light. Versatile.

SMART IMAGING SOLUTIONS

## Bigger Performance. Smaller Package.


The CARESTREAM DRX-1 System unleashed the power of the X-Factor, with the world's first wireless, cassette-sized detector. Its advanced platform allowed a single DRX detector to work seamlessly across the entire family of DRX products. Soon after, we introduced the second-generation DRX detector, a lighter, faster detector with enhanced features.

Now, our newly designed, small-format detector conforms to ISO 4090-size standards and provides greater positioning flexibility, especially in bassinets.

## Advantages At a Glance

The CARESTREAM DRX Plus 2530C Detector offers:

- X-Factor Thinking – enables sharing the detector across our entire DRX family of products.
- A smaller design for fast, easy positioning in pediatric incubator trays.
- Cesium iodide (CsI) design for dose-sensitive pediatric applications.
- Flexible positioning and easy handling, making it ideal for orthopaedic tabletop imaging.
- The ability to share this small-format detector across the entire DRX family of products.
- Smaller pixel pitch for better image quality and improved diagnostic confidence.



# CARESTREAM DRX Plus 2530C Detector



## Additional Detector Benefits:

- Wireless performance, to eliminate the hassles and hazards of cables and reduce the risk of infection.
- Fast performance and instant access to images.
- Compatible with the CARESTREAM DRX-1 Detector Battery Charger to help eliminate additional equipment.
- Adaptable – easily fits into the CARESTREAM DRX-Revolution and CARESTREAM DRX-Revolution Nano Mobile X-ray Systems detector bin.

## Protection Against the Unexpected

The DRX Plus 2530C is eligible for coverage under the CARESTREAM DRX Accident Protection Plan. This uniquely designed plan can be added to your standard Carestream DRX Warranty or Service Agreement. The plan allows replacement of your DRX Plus 2530C Detector at a fraction of the full cost, should damage occur due to an accident in the course of normal clinical use. For you, this means added confidence and peace of mind, convenient, one-call resolution and fast replacement to maximize uptime.

## DRX Plus 2530C Specifications:

Wireless standard: 802.11AC  
Image size: 25 x 30 cm  
Pixel pitch: 0.098 mm  
External dimensions:  
28.2 x 33.2 x 1.5 cm (11.1 x 13.1 x .59 in)



## Industry-Leading Performance: Meet Our Family of DRX Detectors.

Our portfolio of DR detectors includes wireless, shareable and fixed models to meet your demands, based on exam type, detector size, dosage level and budget. Discover the right detector for your needs at [carestream.com/detectors](http://carestream.com/detectors).

[carestream.com/detectors](http://carestream.com/detectors)



"Rx only"

© Carestream Health, Inc., 2019. CARESTREAM is a trademark of Carestream Health. CAT 2000 253 04/19

# Carestream